



# AeroPace® Neurostimulation Therapy System Overview

---



# Difficulty Weaning: Diaphragm Atrophy Occurs Rapidly

Extending time on ventilation and increasing mortality



- ✗ Ventilation weakens diaphragm muscle by more than **50% in less than one day**<sup>1,2</sup>
- ✗ Weakened muscles make it difficult to regain independent breathing, meaning **extended weaning times** and **reduced weaning rates**<sup>1</sup>
- ✗ **Increased likelihood of reintubation**<sup>3</sup> (with a 5-fold increased risk of mortality), **tracheostomy**, and their associated costs<sup>4</sup>
- ✗ **Longer ICU/hospital stays**<sup>3</sup>
- ✗ Leads to increases in respiratory **complications** and **mortality**<sup>3</sup>

1. Dres, M. (2017) *Am J Respir Crit Care Med.*; 195(1): 57-66. doi: [10.1164/rccm.201602-0367OC](https://doi.org/10.1164/rccm.201602-0367OC); 2. Levine, S. (2008) *NEJM*; 358(13): 1327-1355. doi: 10.1056/NEJMoa070447; 3. Zilberberg, M. (2020) *Crit Care Med.*; 48(11): 1587-1594. doi: [10.1097/ccm.0000000000004525](https://doi.org/10.1097/ccm.0000000000004525); 4. Menon, N. (2012) *Respir Care*; 57(10): 1555-1563. doi: [10.4187/respcare.01617](https://doi.org/10.4187/respcare.01617).

# AeroPace vs Standard of Care<sup>1,2</sup>

RESCUE 3 RCT performed across diverse, challenging patient population



UP  
TO **3.2**

**Fewer Days on  
MV in the ICU**



**Reduction in ICU  
Ventilator Time**

UP  
TO **43%**

**More Likely to  
Wean by Day 30**



**Improved Weaning  
Success**

UP  
TO **50%**

**Lower Risk of  
Reintubation\* and  
Tracheostomy**


\*Includes return-to-MV for tracheostomy patients



**Fewer Additional  
Procedures**

UP  
TO **74%**

**Stronger  
Diaphragm**



**Increased  
Strength**




# ➤ **The AeroPace Neurostimulation Therapy System**

- Patient Selection
- Components and Setup
- Neurostimulation Catheter Insertion and Placement
- Mapping and Therapy Sessions
- Super-users: Cleaning, Storage, and MR Compatibility
- Super-users: Troubleshooting (Manual mode and Exclusion)



# AeroPace Therapy: Neurostimulations Delivered through Specialized CVC Exercises and Rehabilitates the Diaphragm

Mechanical ventilation causes diaphragm atrophy, leading to difficulty weaning. AeroPace stimulations contract the patient diaphragm, **strengthening the muscle** and **helping the patient wean** off the ventilator.



# AeroPace® System Safety and Effectiveness Information

## Indications and Contraindications for Use



### Intended Use

The AeroPace System is intended for temporary stimulation of the phrenic nerve(s) to increase diaphragmatic strength.

The AeroPace System is intended for use in hospitals and hospital-type facilities which provide care for patients requiring mechanical ventilation. The device is intended to be used by appropriately trained personnel under the direction of a physician.

### Indication for Use

**The AeroPace System is indicated to improve weaning success – increase weaning, reduce ventilator days, and reduce reintubation – in patients ages 18 years or older on mechanical ventilation  $\geq$  96 hours and who have not weaned.**

### Contraindications

Do not use the AeroPace System with active implanted cardiac pacemakers, defibrillators, or other implantable electronics within proximity to the AeroPace Neurostimulation Catheter. The AeroPace System has not been clinically evaluated for safety with these implantable electronic devices.



- The AeroPace Neurostimulation Therapy System


## ➤ Patient Selection

- AeroPace System Components and Setup
- Neurostimulation Catheter Insertion and Placement
- Mapping and Therapy Sessions
- Super-users: Cleaning, Storage, and MR Compatibility
- Super-users: Troubleshooting (Manual mode and Exclusion)



# AeroPace Therapy Patient Selection

For patients with difficulty weaning from mechanical ventilation







- The AeroPace Neurostimulation Therapy System
- Patient Selection

## ➤ **Components and Setup**

- Neurostimulation Catheter Insertion and Placement
- Mapping and Therapy Sessions
- Super-users: Cleaning, Storage, and MR Compatibility
- Super-users: Troubleshooting (Manual mode and Exclusion)


# The AeroPace System

Three main components integrate system and ventilator



## Neurostimulation Catheter

A temporary, transvenous neurostimulation device with **CVC functionality**, **plus specialized electrodes** that stimulate the phrenic nerves and activate the diaphragm.



## Neurostimulation Console

Controls **adjustable stimulation parameters** that rebuild diaphragm muscle strength and help patients wean from mechanical ventilation.




## Airway Sensor

**Compatible with any ventilator** for single-patient use to deliver AeroPace Therapy with ventilator breaths.




# AeroPace System Setup and Overview



# AeroPace Neurostimulation Catheter Specifications

Includes standard CVC functionality, plus specialized electrodes




**30 neurostimulation electrodes** arranged in two arrays for left and right side phrenic nerve stimulation



- **23cm, 8.5 Fr diameter** fits easily within vein; standard left jugular or left subclavian vein insertion using Seldinger technique
- **Sterile, designed for single-use for up to 30 days.** Follow institutional guidelines for care and maintenance of CVCs.
- When properly placed in central circulation, the Neurostimulation Catheter may be used for **TPN, fluid, and medication delivery.** Infusions may continue during Therapy delivery.
- See Neurostimulation Catheter Kit Instructions for Use for additional information.

High Pressure Lumen (5 ml/sec and flow rate of 42 ml/min) and Guidewire insertion




Medial (yellow) and Proximal (white) Lumen flow rate is 17 ml/min



# AeroPace Neurostimulation Console

Generates stimulations that activate the phrenic nerves



Features of the **Neurostimulation Console** include:

- ✓ **Touchscreen** User Interface
- ✓ **ECG-guided** Catheter insertion
- ✓ Neurostimulation Catheter **Placement** verification
- ✓ **Mapping** to select effective electrodes
- ✓ User-adjustable **Therapy** parameters and notifications for diaphragm activation

# AeroPace Neurostimulation Console

## Buttons and Cables



**Yellow Stop Button:** Pressing this button will quickly stop stimulations. The Console must be powered off after pressing this button.

**Catheter Cable:** The cable with blue connectors electrically connects the AeroPace Neurostimulation Catheter to the Console when in use.

**Airway Sensor Cable:** The cable with dark gray connectors electrically connects the Airway Sensor to the Console.



**Power Button:** Powers the Console on and off.


**Handheld Controller:** The Handheld Controller has a light gray connector. It is used to perform System functions at a distance from the Console.

## Use

- Connect the Airway Sensor to the breathing circuit **before setting up the sterile field** to insert the Neurostimulation Catheter.
- The Airway Sensor connects to the mechanical ventilator at the inspiratory limb of the breathing circuit.
  - *Note: Alternatively, if needed, the Airway Sensor may be placed between the Y-piece and the patient's tracheostomy or endotracheal tube.*
- The **arrow** on the sensor should point toward the patient.
- Is designed to remain connected to the patient's breathing circuit.
- Single-patient use for up to 30 days.

## Cleaning

- Gently **suction** secretions to remove them.
- Water could damage the sensor's electronics; do not rinse it under running water.




# AeroPace Airway Sensor

## Required for System Use in Auto mode



### How it Works

- The Airway Sensor connects to and measures changes in the patient's **breathing circuit pressure**.
- Using **pressure** waveform, the Airway Sensor detects the beginning of ventilator breaths and diaphragm contraction.







- The AeroPace Neurostimulation Therapy System
- Patient Selection
- Components and Setup

## ➤ **Neurostimulation Catheter Insertion and Placement**

- Mapping and Therapy Sessions
- Super-users: Cleaning, Storage, and MR Compatibility
- Super-users: Troubleshooting (Manual mode and Exclusion)

# AeroPace Neurostimulation Catheter Kit

Insert Neurostimulation Catheter before expiration date printed on the kit





- **Neurostimulation Catheter Kit includes:** sterile Neurostimulation Catheter and single-use accessories commonly used to insert a CVC using the Seldinger technique.
- **Does not include:** Ancillary medical supplies such as sterile saline flushes, **local anesthetic**, **skin preparation materials**, sterile gloves, gown, or sterile ultrasound probe cover.



# Insertion Location

- The Neurostimulation Catheter should be inserted into the **left subclavian** or **left jugular** vein.
- Prepare for insertion following the Neurostimulation Catheter Kit **Instructions for Use**, as well as institutional guidelines for the insertion of CVCs using the **Seldinger technique**.



*Precaution: The US Center for Disease Control and Prevention (CDC) recommends a subclavian catheter placement to minimize the risk of infection.*


# Proper Distal Tip Position

## Neurostimulation Catheter tip should be:

- ✓ A few centimeters above the azygous vein or the carina or trachea.
- ✓ In distal third of the SVC and parallel to its wall.
- ✓ About 1 cm above the cavoatrial junction.

**WARNING:** Do not position the Neurostimulation Catheter (or allow it to remain) in the atrium, right ventricle, or too close to the heart. **Incorrect positioning of the Neurostimulation Catheter may result in cardiac arrhythmia or tamponade.**

**Note:** The entire length of the Neurostimulation Catheter from tip to hub is 23 cm. It is **not intended to be inserted all the way to the hub** in most patients.







# Neurostimulation Catheter Insertion

Have a non-sterile assistant available to operate Console

1. Have the non-sterile assistant power on the AeroPace System and select **Catheter Insertion**. You then have the **option to use the ECG** feature as an aid in the Catheter insertion process.
2. If you are using the ECG feature, the sterile clinician will slowly insert the Neurostimulation Catheter using the **Seldinger technique** while watching the screen for signal morphology changes.
  - **Avoid bending, twisting, or torquing the Neurostimulation Catheter** during insertion.
  - **Caution:** Follow all instructions in the *Neurostimulation Catheter Kit Instructions for Use*.




# ECG-Guided Insertion (Optional)

Use clinical assistance to maintain sterile field while operating the Console


## Insert the Neurostimulation Catheter

As the Neurostimulation Catheter first enters the vein over the guidewire, the electrogram will appear flat.




## Look for clear signal while advancing

Advance the Neurostimulation Catheter **slowly**, looking for increased amplitude in the signal. **Stop advancing** upon seeing good definition of the intravascular electrogram. **Retract the guidewire.**





## Do not over-advance (inverted signal)

An inverted/flipped signal morphology indicates that the Neurostimulation Catheter has been **advanced too far** and may be in the atrium; retract the Catheter slowly.



# Cable Clamp

When the Neurostimulation Catheter seems to be in good position, use the **Cable Clamp** to help ensure that the **sterile portion** of the Neurostimulation Catheter remains in place as you prepare to perform the **Placement** procedure using the Console.




# Placement

**After inserting the Catheter, but before suturing it in place, perform Placement on the Console.**

**Placement** will confirm that the Neurostimulation Catheter is in a position that can stimulate the patient's left phrenic nerve.

**How it works:** the System delivers small electrical pulses to electrodes in the Neurostimulation Catheter's **left array** while evaluating data from the Airway Sensor to detect pressure changes that indicate diaphragm contraction.





# Performing Placement

Use assistance to maintain the sterile field

## Placement Instructions




1. Press the **Placement** tab on the Console.
2. Press **Stimulate**.
3. Placement **proceeds automatically**. It typically completes in a few minutes (or less).
4. The System will inform you whether Placement was successful. (If it is unsuccessful, repeat it and/or follow troubleshooting steps as needed.)
5. Once you see the completion screen, follow institutional guidelines to confirm Neurostimulation Catheter tip location (e.g., x-ray).





# Confirming Position of Distal Tip

Locating the Neurostimulation Catheter's **radiopaque marker bands** on an x-ray helps ensure both of the Neurostimulation Catheter's **electrode arrays** are **positioned appropriately** to recruit the patient's phrenic nerves.



**Both Left Array Marker Bands should be left of the mediastinum**

**Both Right Array Marker Bands should be right of the mediastinum**

*Note: X-ray images shown in this presentation have been edited to emphasize Neurostimulation Catheter and Marker Band location. Neurostimulation Catheter location in actual x-rays may be more subtle.*

# Neurostimulation Catheter Too Deep

If the Neurostimulation Catheter tip is too close to the heart, it should be retracted

## Indication the tip is too deep:


- Three marker bands are visible on the right side of the mediastinum.

## Why it matters:

- Neurostimulation Catheter tip may be too close to the heart and cause unwanted cardiac stimulation
- Both electrode arrays might stimulate the right phrenic nerve.

## What to do:

- Neurostimulation Catheter should be retracted approximately 5 cm.



# Neurostimulation Catheter Insertion too Shallow

If the Neurostimulation Catheter is not inserted deep enough, consider replacing it

## Indication the tip is too shallow:


- May be in the brachiocephalic vein.
- Tip angled  $\sim 90^\circ$  to the vertical of the SVC wall: Left sided CVCs, including the Neurostimulation Catheter, should have a tip angle less than  $<40^\circ$  to the vertical in order to avoid abutting the vessel wall.<sup>1</sup>

## Why it matters:

- May not be able to stimulate the right phrenic nerve.
- Risk of vessel wall damage or perforation.

## What do to:

- Consider replacing the Neurostimulation Catheter over the wire, placing the tip deeper.



1. Stonelake, P. (2006) *British Journal of Anaesthesia*; 96 (3): 335–40.

# Suturing: 20-23cm Depth

**20-23cm**

Suture hub

( $\leq 20\text{cm}$  hub+clamp)

Sutures

Fully secure the Neurostimulation Catheter (with or without the secondary clamp, depending on its depth) taking care to **minimize acute bends in the Neurostimulation Catheter** where it meets the hub.




*Note: Images shown are of simulated bench model.*

*Note: Always follow institutional guidelines for infection prevention and dressing of the Neurostimulation Catheter insertion site.*




# Suturing: $\leq 20\text{cm}$ Depth



Assembled  
Secondary  
Clamp

1. **For insertions  $\leq 20\text{cm}$  only**, fit the blue Neurostimulation Catheter clamp onto the shaft of the Neurostimulation Catheter as close to the Neurostimulation Catheter's insertion point as possible.
2. Press the white fastener over the blue clamp.
3. Suture the clamp assembly to the patient's skin using its holes. *Never directly suture into the Neurostimulation Catheter shaft.*
4. Suture the blue Neurostimulation Catheter hub to the patient's skin using its holes.

Insertion Point




*Note: Always follow your institutional guidelines for infection prevention and dressing of the Neurostimulation Catheter insertion site.*



# Care After Insertion

- The Neurostimulation Catheter is designed for **single insertion and use**.
- Follow institutional guidelines for care and maintenance of CVCs.
- The Neurostimulation Catheter incorporates electrical neurostimulation capabilities and must be handled with care. **Avoid bending, twisting, or torquing the Neurostimulation Catheter**, including during dressing changes.





- The AeroPace Neurostimulation Therapy System
- Patient Selection
- Components and Setup
- Neurostimulation Catheter Insertion and Placement

## ➤ **Mapping and Therapy Sessions**

- Super-users: Cleaning, Storage, and MR Compatibility
- Super-users: Troubleshooting (Manual mode and Exclusion)


# AeroPace Therapy

60 stimulations twice daily



## AeroPace Therapy Overview

- Power on the Console and Select **Therapy Session**.
- **The Console will prompt you to perform Mapping** before of each Therapy session, to select suitable electrodes.
- One Therapy session is **60 stimulations**, titrated up to the **maximum intensity** possible while ensuring patient comfort.
- Therapy is provided **twice daily**.




# Mapping

Selects electrodes that effectively contract the diaphragm

**Before each Therapy session, you will perform Mapping.**

**Mapping** identifies and selects Neurostimulation Catheter electrodes that are positioned to effectively contract the diaphragm on both sides, and determines the minimum electrical charge required for diaphragm contraction.





**Tiny pressure deflections** or “twitches” detected by the Airway Sensor during Mapping allow the System to **identify suitable electrodes** for Therapy, and determine the **baseline Therapy Levels** needed.





# Mapping Screen

Selects electrodes that effectively contract the diaphragm



**Enable/Disable**  
Enables or disables entire screen

**Auto/Manual**  
(See Troubleshooting section)

**Exclusion**  
(See Troubleshooting section)

**End Session**  
Return to main menu

**Notification Box**  
Displays contextual System and Mapping process information

**Stop Button**  
Stops the delivery of stimulations if needed

**Stimulate Button**  
Press to begin Mapping

**Hemidiaphragm Selectors**  
Selects which side(s) to Map. The default setting is to Map both sides

**Progress Bars**  
Shows the progression of Mapping process on each side

The screenshot shows the 'Mapping' screen of the Lungpacer device. At the top, there is a status bar with 'B.X.XXXX.XX.X.XXXX.XX' and a timestamp '19-Apr-2020 16:00'. Below this is a control bar with a green 'Enabled' toggle, 'AUTO' and 'MANUAL' buttons, and tabs for 'Mapping', 'Therapy', 'Exclusion', and 'End Session' (with an 'X' icon). The main area features a central illustration of lungs with 'Right' and 'Left' labels. On either side of the lungs are 'Right Side' and 'Left Side' selectors, each with a checked checkbox and a progress bar below it. A notification box at the bottom of the main area contains the text 'Both sides selected. Press Stimulate to begin Mapping.' At the very bottom are a large blue 'Stimulate' button with a play icon and a grey 'Stop' button with a square icon.

# Performing Mapping



Mapping proceeds and completes automatically

## Mapping Instructions

1. After Placement completes successfully, the **Mapping** tab automatically displays.

*Note: If performing Mapping after powering on the Console, select **Therapy Session** from the main menu to navigate directly to the Mapping tab.*

2. If anatomical restrictions limit Therapy to only one side, deselect the relevant side.
3. Press **Stimulate** to begin.
4. **Mapping completes automatically** on the selected side(s). It typically completes in a few minutes (or less). **The System will inform you whether Mapping was successful.** (If it is unsuccessful, repeat it and/or follow troubleshooting steps as needed.)



# Therapy Overview

Once Mapping is complete, you are ready to begin **Therapy**.

The goal of Therapy is to deliver **60 stimulations** at the **maximum Therapy Levels** possible while still ensuring patient comfort.



# Therapy Screen: Stimulation Controls

## Stim Frequency

The rate of electrical pulses within a stimulation. ***This should not be changed***

## Skip Breaths

The number of ventilator breaths between each stimulation

## Stim Duration

Set as long as possible to match ventilator Inspiratory Time

## Stims Programmed

The number of stimulations (0-10) delivered with each press of the Stimulate button

## Therapy Levels (Left and Right)

Increase (+) or decrease (-) as needed on each side to optimize Therapy

## Total Stims

The number of stimulations delivered during the Therapy session so far

## Settings


Access a pop-up where Therapy session settings can be modified


## Stop Button

Stops the delivery of stimulations if needed

## Stimulate Button


Press to deliver the set number of Stims Programmed







To begin and titrate Therapy:

1. Set **Stim Duration** as long as possible while matching patient's **Inspiration Time** on ventilator.
2. Begin to titrate Therapy by pressing **Stimulate** to deliver initial stimulations.
3. Observe patient, and **increase the Therapy Level** on both the left and right side if possible, by pressing the **plus (+) buttons on each side several times**.
4. Press **Stimulate** to deliver another stimulation (or two) and observe the patient.
5. Continue this process until reaching the maximum possible Therapy Level on each side that can be delivered to the patient comfortably.







To complete Therapy:

1. Continue to deliver stimulations until the **Total Stims** counter in the upper right-hand corner of the screen reaches 60.

*Note: You may press Settings and increase **Stims Programmed** to any value from 1 to 10 as you approach 60 stimulations.*

2. Once **60 stimulations** are delivered, press **End Session** to display session data.
- **Power off the System, and disconnect the Catheter and Airway Sensor** until the next Therapy session. The cables may remain attached to the Console, neatly arranged on the handle.


**Large pressure deflections** during breaths indicate AeroPace stimulations are contracting the patient diaphragm to exercise it, and generating negative pressure.



**Incrementally increase the Therapy Levels on each side to reach the maximum possible while maintaining patient comfort, to deliver 60 stimulations.**

# AeroPace Therapy Patient Selection

For patients with difficulty weaning from mechanical ventilation





# Knowledge Check

## AeroPace System

Scan the QR code with your phone to begin.





- The AeroPace Neurostimulation Therapy System
- Patient Selection
- Components and Setup
- Neurostimulation Catheter Insertion and Placement
- Mapping and Therapy Sessions


## ➤ **Super-users: Cleaning, Storage, and MR Compatibility**

- Super-users: Troubleshooting (Manual mode and Exclusion)



# Cleaning

Clean the Console, Airway Sensor Cable, Catheter Cable, and Handheld Controller using Virox, Cidex, or CaviWipes, following your institutional policy for cleaning bedside equipment.



## To store the AeroPace System:

- Disconnect cables
- Ensure the System is stored within the following environmental conditions:
  - Console and Cables Temperature Range: -10 to 60 °C (14 to 140 °F)
  - Console and Cables Relative Humidity Range: 15% to 90%, non-condensing
  - Console and Cables System Altitude: Up to 8000 ft (2438 m) above sea level
  - Neurostimulation Catheter Temperature Range: 18 to 28 °C (64.5 to 82.4 °F)




## Zone 1

Top of head to cervical region, C4-5

**Catheter Cable must be oriented toward feet**

SAR Limits:  
WBA SAR, 0.1 W/kg  
Head SAR, 0.2 W/kg



## Zone 2

Cervical region, C4-5  
to below knee

## Zone 3

Below knee to  
bottom of foot

## Zone 1


Top of head to cervical region, C4-5

## Zone 2

Cervical region, C4-5  
to below knee

**Catheter Cable  
must be oriented  
toward head**

SAR Limits:  
WBA SAR, 0.1 W/kg



## Zone 3

Below knee to bottom of foot

**Catheter Cable must be  
oriented toward head**

SAR Limits:  
WBA SAR, 2.0 W/kg  
(Normal Operating Mode)

- **The AeroPace Neurostimulation Catheter is MR Conditional.** A patient may be safely scanned under conditions listed in the Instructions for Use and with the Neurostimulation Catheter oriented as shown.
- **Warning:** The AeroPace Neurostimulation Catheter is MR Conditional, but the Neurostimulation Console and Cables are MR Unsafe. Always disconnect the Catheter from the Neurostimulation Console before conducting MRI, and keep the Console away from MR equipment.



- The AeroPace Neurostimulation Therapy System
- Patient Selection
- Components and Setup
- Neurostimulation Catheter Insertion and Placement
- Mapping and Therapy Sessions
- Super-users: Cleaning, Storage, and MR Compatibility

## ➤ **Super-users: Troubleshooting (Manual mode and Exclusion)**

Placement and Mapping


Auto Mode	Manual Mode
Diaphragm contractions detected by <b>Airway Sensor</b>	Diaphragm contractions detected by <b>clinician palpating patient</b> (and then providing feedback to Console).
Stimulations delivered when <b>Airway Sensor</b> detects patient breaths	Stimulations delivered at end-expiration at set <b>Stim Rate</b> , which must be set to match ventilator breath rate to avoid dyssynchrony.

Therapy




## Placement Troubleshooting

- ☐ If Placement is unsuccessful, repeat Placement in Auto mode.
- ☐ Consider repositioning the patient.
- ☐ Assess the Neurostimulation Catheter's location and consider repositioning it.
- ☐ If unwanted effects occur, use the Stop button. Use the electrode Exclusion function if unwanted effects continue.
- ☐ Try Manual mode if Auto mode is persistently unsuccessful or results in unwanted effects.



# Manual Placement

To perform manual Placement, press Stimulate while you **palpate the patient to feel for diaphragm contractions**, and then provide feedback to Console using the on-screen buttons or Handheld Controller.




Select **Unsure** if you're not sure if you felt a contraction or not.

Press **Stimulate** if you did not feel anything.

Select **Confirm** if you felt a diaphragm contraction.


## Mapping Troubleshooting

- ☐ If Auto Mapping is unsuccessful, repeat it.
- ☐ Consider repositioning the patient.
- ☐ Try Manual mode if Auto Mapping is persistently unsuccessful. See the System User Guide for details.
- ☐ If unwanted effects occur, use the Stop button. Use the electrode Exclusion function if unwanted effects continue.



# Manual Mapping

To perform manual Mapping, press Stimulate while you **palpate the patient to feel for diaphragm contractions**, and then provide feedback to Console using the on-screen buttons or Handheld Controller.




Select **Unsure** if you're not sure if you felt a contraction or not.

Press **Stimulate** if you did not feel anything.

Select **Confirm** if you felt a diaphragm contraction.

## Therapy Troubleshooting

- ❑ If patient discomfort, muscle twitching, or other unwanted effects occur:
  - ❑ Use the **Stop** button to quickly stop stimulations if needed.
  - ❑ Incrementally decrease the **Therapy Level** on affected side.
  - ❑ Use electrode **Exclusion** only if needed.
- ❑ If a persistent **Double Trigger Warning** occurs: Increase number of Skip Breaths and/or reduce the Stim Duration if clinically possible.






During manual Therapy, Stimulations are delivered at the set Stim Rate, which must be configured to match the ventilator breath rate. Additionally, the Stimulate button must be pressed at end-expiration of the patient's breathing cycle to avoid ventilator dyssynchrony.

## To provide Manual Therapy:


1. Adjust the Therapy Levels to optimize Therapy.
2. Set the **Stim Rate** and **Stim Duration** carefully to match the patient's ventilator settings, and press **Stimulate at end-expiration**,
3. Continue to assess the patient during Therapy.
4. Deliver 60 total stimulations (2x/day).



# Exclusion

Use the Exclusion tab to identify and then temporarily disable (exclude) the specific Catheter electrodes that cause unwanted effects. Exclusion can be performed as part of Placement, Mapping, or Therapy.

Press Stimulate while you **observe the patient to watch for unwanted effects**, and then provide feedback to Console using the on-screen buttons or Handheld Controller.



Select **Unsure** if you're not sure if you saw an unwanted effect.

Press **Stimulate** if you did not observe anything of note.

Select **Confirm** if you observed an unwanted effect.



260 Sierra Drive, #116  
Exton, PA 19341

**[Lungpacer.com](https://Lungpacer.com)**