



Handheld Controller

Lungpacer Handheld Controller

A component of the Lungpacer Diaphragm Pacing Therapy System

Instructions for Use

Rx Only

For use only with the Lungpacer Diaphragm Pacing Therapy System (DPTS) under an Emergency Use Authorization. Refer to the Lungpacer Diaphragm Pacing Therapy System (DPTS) Instructions for Use for information.



Lungpacer Medical, Inc.
130 – 601 W Cordova Street
Vancouver, BC V6B 1G1
Canada
+1 (778) 655-2100
www.lungpacer.com

REF 200-0029


IFU Artwork: 650-0179 Rev A
Effective Date: 2020-04-17



1. DEVICE DESCRIPTION

The Handheld Controller is a component of the Lungpacer Diaphragm Pacing Therapy System (DPTS). The Handheld Controller connects to the Lungpacer Control Unit.

A diagram of the Handheld Controller is shown below:



Handheld Controller

The Lungpacer DPTS is intended for use by qualified, trained personnel under the direction of a physician.

The Handheld Controller connects to the Lungpacer Control Unit (LCU) to allow remote stimulation of the left and/or right phrenic nerves causing the diaphragm to contract. The Connector at the end of the Handheld Controller is intended only for connection to the LCU.

2. INTENDED USE

Refer to Lungpacer Diaphragm Pacing Therapy System Instructions for Use.

3. DISINFECTION INSTRUCTIONS

Refer to Lungpacer Diaphragm Pacing Therapy System Instructions for Use.

4. SYMBOL GLOSSARY



Manufacturer



Catalogue number



Caution



Consult instructions for use



Keep dry



Keep away from sunlight



Do not use if package is damaged



Prescription only