



Intermediate Cable

Lungpacer Intermediate Cable

A component of the Lungpacer Diaphragm Pacing Therapy System.

Instructions for Use

Rx Only

For use only with the Lungpacer Diaphragm Pacing Therapy System (DPTS) under an Emergency Use Authorization. Refer to the Lungpacer Diaphragm Pacing Therapy System (DPTS) Instructions for Use for information.



Lungpacer Medical, Inc.
130 – 601 W Cordova Street
Vancouver, BC V6B 1G1
Canada
+1 (778) 655-2100
www.lungpacer.com

REF 200-0031


IFU Artwork: 650-0180 Rev A
Effective Date: 2020-04-17



1. DEVICE DESCRIPTION

The Intermediate Cable is a component of the Lungpacer Diaphragm Pacing Therapy System (DPTS). The Intermediate Cable connects between the Lungpacer Control Unit and a Lungpacer LIVE Catheter.

A diagram of the Intermediate Cable is shown below:



Intermediate Cable

The Lungpacer DPTS is intended for use by qualified, trained personnel under the direction of a physician.

The Intermediate Cable connects the Lungpacer Control Unit (LCU) to the LIVE Catheter. The larger connector with the Lungpacer logo at one end of the Intermediate Cable is intended only for connection to the LCU. The smaller connector with the patient icon at the other end of the Intermediate Cable is intended only for connection to a LIVE Catheter.

2. INTENDED USE

Refer to the Lungpacer Diaphragm Pacing Therapy System (DPTS) Instructions for Use.

3. DISINFECTION INSTRUCTIONS

Refer to the Lungpacer Diaphragm Pacing Therapy System (DPTS) Instructions for Use.

4. SYMBOL GLOSSARY



Manufacturer



Catalogue number



Caution



Consult instructions for use



Keep dry



Keep away from sunlight



Do not use if package is damaged



Prescription only